




**Amerisoft Medical Products**  
*Quality Patient Care Management Products*

## Deluxe Vinyl Padded Half Rail Bed Safety Rails Pads



Deluxe vinyl padded bed safety rail pads help prevent patient injury from bed safety rails. Manufactured with high-density 1" foam (CA flame retardant) provides soft protection, with heavy-duty blue vinyl cover which easily wipes clean. Each pad includes two heavy duty nylon safety straps with deluxe heavy-duty acetal buckles to keep pad securely in place. Fits most standard half length and split rail length long-term care and home care style bed rails. Measures 29" long by 14" high. Sold by the pair or in a pack of four (4).

**Reorder  
Number:**

**Description:**

101029

29" Deluxe vinyl bed rail pads. Sold by the pair.

101029-4

29" Deluxe vinyl bed rail pads. Four (4) per package. Sold by the pack.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.