



Amerisoft Medical Products
Quality Patient Care Management Products

Deluxe Vinyl Padded Bed Safety Rails Pads



Deluxe vinyl padded bed safety rail pads help prevent patient injury from bed safety rails. Manufactured with high-density 1" foam (CA flame retardant) provides soft protection, with heavy-duty blue vinyl cover which easily wipes clean. Each pad includes four heavy duty nylon safety straps with deluxe heavy-duty acetal buckles to keep pad securely in place. Fits most full length long-term care and home care style bed rails. Available in 60" and 72" overall lengths. Packaged individually. Sold by the pair.

**Reorder
Number:**

Description:

101060

60" Deluxe vinyl bed rail pads. Sold by the pair.

101072

72" Deluxe vinyl bed rail pads. Sold by the pair.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.