



Amerisoft Medical Products
Quality Patient Care Management Products

Deluxe Quilted Padded Half Rail Bed Safety Rail Pads



Deluxe quilted bed safety rail pads help prevent patient injury from bed safety rails. Manufactured with pillowy soft polyfiber fill and an attractive navy blue rip-stop nylon cover that is both water-resistant and bacteriostatic. Provides soft, yet effective protection. Each rail pad features two heavy-duty nylon straps with high-strength, dual quick-release plastic buckles. Fits most standard long term care and home care style 1/2 bed rails. Meets or exceeds flame retardancy requirements of California Technical Bulletins 116 & 117. Available in pairs or sets of four (4) pads.

Reorder

Number:

Description:

101229

Deluxe quilted padded 1/2 bed rail pads. Sold by the pair.

101229-4

Deluxe quilted padded 1/2 bed rail pads. Sold by the set of four (4) pads.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.