



Amerisoft Medical Products
Quality Patient Care Management Products

Valore Line Padded Bed Safety Rail Pads

NEW!



Valore line padded full length bed rail pads. Manufactured with a 1.5" foam core and durable light blue bacteriostatic vinyl cover that is both waterproof and flame retardant. Includes four (4) heavy-duty vinyl straps with hook & loop closure to secure pads to bed rails. Full length zipper allows for easy removal of cover. Cover may be cleaned with liquid disinfectant. Measures 1.5" x 13" x 72". Sold by the pair.

Reorder

Number:

101485

Description:

Valore line vinyl bed rail pads. 1.5" x 14" x 72" Sold by the pair.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.