



Amerisoft Medical Products
Quality Patient Care Management Products

Full Leg Abduction Wedge Cushion



Helps to prevent leg abduction and contracture for bed-bound patients, while keeping the patient as comfortable as possible. Manufactured with a medium density foam core, surrounded by soft, breathable polyfiber fill. Covered with a soft, comfortable, breathable baby blue velour fleece material. Secured with two soft fabric straps with hook and loop fasteners for extra stability. Full-length nylon strap with four loops helps to keep hook and loop retaining straps in the desired location. Measures 22" long, by 7" wide at the knee to 14" wide at the ankle. Zippered closure allows for removal of cover for washing. Cover is machine washable at low temperatures. Sold by the each.

**Reorder
Number:**

335000

Description:

Full Leg Abduction Wedge Cushion. Sold by the each.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.