



**Amerisoft Medical Products**  
*Quality Patient Care Management Products*

## **Foot Drop & Blanket Support Wedge for Bed**



Foot drop wedge support for bed. Helps prevent foot drop, while keeping sheets and bed linens off resident's feet. Manufactured with high density foam fill and heavy-duty, reinforced, bacteriostatic gray vinyl that easily wipes clean with soap & water or liquid disinfectant. Can use either straight or angled surface for foot positioning needs. Cushion is held securely in position with two extra-long, high strength nylon securing straps - excess length can be cut once support is secured to bed (extra long strap kits for full, queen and king size mattresses are available separately). Quick release buckles are provided on both sides of cushion to accommodate adjustment from either side of the bed. Built-in handle affords easy transport. Measures 13" high by 24" wide, with a gentle slope from 10" thick on the bottom to 7" thick on the top of the cushion.

**Reorder  
Number:**

401300

**Description:**

Foot Drop & Blanket Support Wedge for Bed. Sold by the each.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.