



Amerisoft Medical Products
Quality Patient Care Management Products

Foam Wheelchair Cushions



Medium density foam fill provides durable support and excellent pressure reduction at an economical price. Improved water-resistant dark blue nylon cover with waterproof ski-resistant backing. Meets or exceeds flammability requirements of California bulletins 116 and 117. Available in a variety of sizes. Sold by the each.

**Reorder
Number:**

Description:

421003	18" x 16" x 3". Sold by the each.
421004	18" x 16" x 4". Sold by the each.
421013	20" x 18" x 3". Sold by the each.
421014	20" x 18" x 4". Sold by the each.
421023	22" x 18" x 3". Sold by the each.
421024	22" x 18" x 4". Sold by the each.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.