



Amerisoft Medical Products
Quality Patient Care Management Products

Valore Line Wheelchair Foam Wedge Cushion



Created to offer better value for the budget-minded, Amerisoft's new Valore line offers durable performance at a value price. The new Valore foam wheelchair wedge cushion provides pressure reduction and a classic wedge shape, rising from 1.5" in the rear of the cushion to 3.5" in the front of the cushion to help prevent forward sliding. Manufactured with a soft foam core and a durable, incontinent-resistant, bacteriostatic gray vinyl cover. Meets or exceeds flammability requirements of California Technical Bulletins 116 and 117. Measures 18" x 16" x 1.5" to 3.5". Integrated hook & loop straps help keep cushion firmly in place. Custom colors and sizes are available - contact your local distributor or our sales department for more information. Sold by the each.

**Reorder
Number:**

421090

Description:

Valore Line Wedge cushion. 18" x 16" x 1.5" to 3.5". Sold by the each.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.