



## **Padded Wheelchair Comfort Armrest Cushions**



Wheelchair armrest cushions provide padded comfort, protection and pressure relief for the arms and elbows of wheelchair users. Manufactured from a 1" high density foam fill (CA flame retardant) for optimum pressure relief, and covered with a soft, comfortable, durable black ripstop nylon. Includes full length hook and loop fasteners to help keep cushion securely in place. Sold by the pair.

### **Reorder**

### **Number:**

### **Description:**

- |        |  |
|--------|--|
| 851010 | Half length padded wheelchair armrest cover. 11". Dark blue. Sold by the pair. |
| 851011 | Half length padded wheelchair armrest cover. 11". Black. Sold by the pair.     |
| 851012 | Full length padded wheelchair armrest cover. 16". Dark blue. Sold by the pair. |
| 851013 | Full length padded wheelchair armrest cover. 16". Black. Sold by the pair.     |

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.