



Amerisoft Medical Products
Quality Patient Care Management Products

Lateral Side Supports for Wheelchair & Gerichair



Lateral side supports provide comfortable, snug support and are effective for promoting proper sitting posture. Manufactured from durable fiberboard with high density foam fill (CA flame retardant) covered with a heavy-duty, wipe-clean blue vinyl cover for durability. Single support helps to control lateral leaning, while use of two supports helps to control bilateral leaning. Includes hook & loop closures to help keep support securely in place. Measures 21" long x 16" high. Sold by the pair.

**Reorder
Number:**

852010

Description:

Lateral side supports for wheelchair & gerichair. Sold by the pair.

Need custom sizes or colors? We can help! Contact your authorized Amerisoft Medical distributor, or our Customer Service Dept. for more info.